Twisted-files Adaptive: Un innovador enfoque en la instrumentación Endodóntica NiTi

Dr. Adi Moran

Twisted-files Adaptive: Un innovador enfoque en la instrumentación Endodóntica NiTi

Dr. Adi Moran

DMD

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References

1. Hummonen S, Orstavik D. Radiological aspect of apical periodontitis. *Endod Topics* 2002;1:3-25
2. Yared G. Canal preparation using only one Ni-Ti rotary instrument: preliminary observations. *Int Endod J* 2008;41:339-44.
3. phase-transformation temperatures of TF were significantly higher than those of machined files, with the bending load values being significantly lower for these files in the elastic and SE ranges (Hou et al. 2011). Hou XM, Yahata Y, Hayashi Y, Ebihara A, Hanbawa T, Suda H (2011) Phase transformation behavior and bending property of twisted nickel-titanium instruments. *International Endodontic Journal* 44, 253-8.
4. Hulsmann M, Peters OA, Dummer PMH (2005) Mechanical preparation of root canals, shaping goals, techniques and means. *Endodontic Topics* 10, 30-76.
5. Bahia MGA, Bueno VCL (2005) Decrease in the fatigue resistance of nickel-titanium rotary instruments after clinical use in curved canals. *Oral Surgery Oral Medicine Oral Pathology Oral Radiology and Endodontology* 100, 249-55.
6. Anderson ME, Price JW, Parashos P (2007) Fracture resistance of electropolished rotary nickel-titanium endodontic instruments. *Journal of Endodontics* 33, 1212-6.
7. Gambarini G, Gerosa R, De Luca M, Garala M, Testarelli L (2008a) Mechanical properties of a new and improved nickel-titanium alloy for endodontic use: an evaluation of file flexibility. *Oral Surgery Oral Medicine Oral Pathology Oral Radiology and Endodontology* 105, 798-800.
8. Otsuka K, Ren X (2005) Physical metallurgy of Ti-Ni-based shape memory alloy. *Progress in Materials Science* 50, 511- 678.
9. Burklein, Dr med dent, and Edgar Schafer, Prof Dr med dent. Apically Extruded Debris with Reciprocating Single-File and Full-sequence Rotary Instrumentation Systems Sebastian JOE — Volume 38, Number 6, June 2012 850-852
10. Gambarini1, R. Gergi2, A. Naaman2, N. Osta2 & D. Al Sudani3. Cyclic fatigue analysis of twisted file rotary NiTi instruments used in reciprocating motion. *Int Endod Journal* 2012 1-5.
11. Vadhana S, SaravanaKarthikeyan B, Nandini S, Velmurugan N. Cyclic fatigue resistance of RaCe and Mtwo rotary files in continuous rotation and reciprocating motion. *J Endod* 2014;40:995-9. 14.
12. Gavini G, Caldeira CL, Akisue E, et al. Resistance to flexural fatigue of Reciproc R25 files under continuous rotation and reciprocating movement. *J Endod* 2012;38: 684-7. 15.
13. Goldberg M, Dahan S, Machtou P. Centering ability and influence of experience when using WaveOne single-file technique in simulated canals. *Int J Dent* 2012; 2012:206321.
14. Lim YJ, Park SJ, Kim HC, Min KS. Comparison of the centering ability of Wave.One and Reciproc nickel-titanium instruments in simulated curved canals. *Restor Dent Endod* 2013;38:21-5
15. Marzouk AM, Ghoneim AG. Computed tomographic evaluation of canal shape instrumented by different kinematics rotary nickel-titanium systems. *J Endod* 2013; 39:906-9.
16. Junaid A, Freire LG, da Silveira Bueno CE, et al. Influence of single-file endodontics on apical transportation in curved root canals: an ex vivo micro-computed tomographic study. *J Endod* 2014;40:717-20.
17. Wan J, Rasimick BJ, Musikant BL, Deutsch AS. A comparison of cyclic fatigue resistance in reciprocating and rotary nickel-titanium instruments. *Aust Endod J* 2011; 37:122-7.
18. You SY, Bae KS, Baek SH, et al. Lifespan of one nickel-titanium rotary file with reciprocating motion in curved root canals. *J Endod* 2010;36:1991-4
19. Burklein S, Schafer E. Apically extruded debris with reciprocating single-file and full-sequence rotary instrumentation systems. *J Endod* 2012;38:850-2. 91.
20. Burklein S, Hinschitzka K, Dammaschke T, Schafer E. Shaping ability and cleaning effectiveness of two single-file systems in severely curved root canals of extracted teeth: Reciproc and WaveOne versus Mtwo and ProTaper. *Int Endod J* 2012;45: 449-61.
21. Çapar ID, Arslan H, Akcay M, Ertas H An in vitro comparison of apically extruded debris and instrumentation times with ProTaper Universal, ProTaper Next, Twisted File Adaptive, and HyFlex instruments. *Journal of Endodontics* 2014 40, 1638-41.
22. Gambarini, D. AL Sudani, S. DI Carlo, G. Pompa, A. Pacifici. I. Pacifici. L. Testarelli. Incidence and intensity of postoperative pain and periapical inflammation after Endodontic treatment with two different instrumentation techniques. *Euro J. Inflam* 2012; 10:99-103
23. Gambarini G, Gergi R, Naaman A, et al. Cyclic fatigue analysis of twisted file rotary NiTi instruments used in reciprocating motion. *Int Endod J* 2012;45:802-6.
24. Thompson SA. An overview of nickel-titanium alloys used in dentistry. *International Endodontic Journal*, 33, 297-310, 2000.
25. Perez-Higueras JJ, Arias A, de la Macorra JC. Cyclic fatigue resistance of K3, K3XF, and twisted file nickel-titanium files under continuous rotation or reciprocating motion. *J Endod* 2013;39:1585-8.
26. Ordinola-Zapata R, Bramante CM, Duarte MA, et al. Shaping ability of reciproc and TF adaptive systems in severely curved canals of rapid microCT-based prototyping molar replicas. *J Appl Oral Sci* 2014;22:509-15.
27. Hamid Jafarzadeh, Paul V. Abbott. Ledge Formation: Review of a Great Challenge in Endodontics. *Journal of Endodontics*, Volume 33, Issue 10, October 2007, Pages 1155-1162
28. Gutmann JL, Gao Y. Alteration in the inherent metallic and surface properties of nickel-titanium root canal instruments to enhance performance, durability and safety: a focused review. *International Endodontic Journal*, 45, 113-128, 2012.
29. De-Deus G, Barino B, Zamolyi RQ, et al. Suboptimal debridement quality produced by the single-file F2 ProTaper technique in oval-shaped canals. *J Endod* 2010;36: 1897-900.
30. Dietrich MA, Kirkpatrick TC, Yaccino JM. In vitro canal and isthmus debris removal of the self-adjusting file, K3, and WaveOne files in the mesial root of human mandibular molars. *J Endod* 2012;38:1140-4.
31. H Kim, J Jum, B Hur, G S Cheung. Cyclic fatigue and fracture characteristics of ground and twisted Nickel-Titanium rotary files. *J Endod* 2010; 36: 147-152.
32. R Gergi, R Arbab-Chirani, N Osta, A Naaman. Micro-Computed tomographic evaluation of canal transportation instrumented by different kinematics rotary nickel-titanium instruments. *Int Endod J* 2015; 48: 162-170.
33. A Dhingra, Neetika. Glide path in Endodontics. *Endodontology* 2014;Vol 26, 1, 217-222.
34. Berutti E1, Paolino DS, Chiandussi G, Alovisi M, Cantatore G, Castellucci A, Pasqualini D. Root canal anatomy preservation of WaveOne reciprocating files with or without glide path. *J Endod.* 2012 Jan;38(1):101-4.