A background pattern of light gray hexagons arranged in a honeycomb-like structure, with some hexagons missing to create a grid-like effect. A solid black horizontal bar is positioned across the middle of the page, containing the product name and description.

Vertise™ *Flow*

Self-Adhering Composite

Directions For Use

Vertise™ Flow

Directions For Use

Vertise™ Flow composite is a self-adhesive, light-cured, resin-based composite dental restorative designed for direct placement. By incorporating Kerr's OptiBond adhesion technology, Vertise™ Flow composite eliminates the additional steps of etching/priming/bonding otherwise necessary to bond a resin composite to dentin and enamel. Vertise™ Flow composite offers high bond strength, high mechanical strength, and other physical attributes comparable to traditional flowable composites. Vertise™ Flow composite is indicated for small Class I and base/liner for Class I and II restorations. Additional indications include pit and fissure sealant, repair of enamel defects, repair of porcelain restoration, blocking of undercut, minor occlusal build-ups in non-stress bearing areas, and incisal abrasions.

Vertise™ Flow composite is available in a variety of shades. Shade is indicated on the label of each syringe.

Storage:

- Store between 2-8°C (36 -46°F).
- Refrigerate immediately upon receipt and when not in use.
- Allow material to reach room temperature prior to use.
- Expiration date is marked on all packaging. Do not use after expiration date.

To stop unwanted paste flow from tip, pull back the plunger after dispensing.

It is recommended that a reliable, high-quality, visible light source (such as Kerr's Demi, L.E.Demetron II, or Optilux 501) be used to cure this product.

Shield Vertise™ Flow composite from visible light prior to use.

PRIOR TO PLACEMENT

Recommendations on proper placement

- Isolate tooth from adjacent teeth. Use clear thin strips. Rubber dam is highly recommended.
- Dry dentin surfaces are shown to enhance bond strengths.
- Vertise™ Flow composite is self-adhesive and therefore does not require etching or bonding protocol prior to placement.
- Placement sequences, such as those outlined in this instruction booklet, should be followed closely.
- Take a moment to identify the surface to which you are bonding.
- Please take care to ensure that your air line is free of oil and other contaminants.

PLACEMENT OF VERTISE™ FLOW COMPOSITE

For small Class I restorations

1. Prepare Cavity. Wash thoroughly with water spray and air dry at maximum air pressure for 5 seconds. Note: Whenever possible, prep (bevel) all margins prior to application of Vertise™ Flow.
2. Select the desired shade.
3. Dispense Vertise™ Flow onto prep with provided dispensing tip. Use provided brush to apply Vertise Flow to the entire cavity wall and beveled area with moderate pressure for 15-20 seconds to obtain a thin layer (<0.5mm). Remove excess material around margins with the brush if necessary.
Note: Replace syringe cap after each use to prevent the resin from hardening in the syringe.
4. Light cure for 20 seconds†. For A3.5 and Universal Opaque shades, light cure each increment for 40 seconds.

5. After lining the cavity wall and beveled area with Vertise™ Flow, build the restoration with more Vertise™ Flow in increments of 2mm or less.
6. Light cure each increment for 20 seconds†. For A3.5 and Universal Opaque shades, light cure each increment for 40 seconds.
7. Finish and polish using a polishing system such as KerrHawe OptiDisc, Occlubrush or the HiLuster.

For Class I and Class II restorations as Total Prep Liner

1. Prepare Cavity. Wash thoroughly with water spray and air dry at maximum air pressure for 5 seconds. Note: Whenever possible, prep (bevel) all margins prior to application of Vertise™ Flow.
2. Select the desired shade.
3. Dispense Vertise™ Flow onto prep with provided dispensing tip. Use provided brush to apply Vertise Flow to the entire cavity wall and beveled area with moderate pressure for 15-20 seconds to obtain a thin layer (<0.5mm). Remove excess material around margins with the brush if necessary.
4. Light cure for 20 seconds†. For A3.5 and Universal Opaque shades, light cure each increment for 40 seconds.
5. After lining the cavity wall and beveled area with Vertise™ Flow, build the restoration with a traditional universal composite resin such as Premise™ or Herculite™ Ultra. The thickness of each increment of the traditional universal composite should be no more than 2mm.
6. Light cure each increment according to manufacturer's recommendation.
7. Finish and polish using a polishing system such as KerrHawe OptiDisc, Occlubrush or the HiLuster.

For Pit and Fissure Sealant

1. Isolate the area with a rubber dam or cotton rolls.
2. Prior to sealant application, clean occlusal surfaces as follows:
 - a) Air-abrade. Rinse thoroughly and dry.OR
 - b) Clean with a non-fluorinated pumice first, then apply a 37.5% phosphoric gel etch to the tooth for 15 seconds. Rinse for 10 seconds. Dry for 5 seconds. Etching is necessary for cleaning purposes only.
3. Dispense Vertise™ Flow onto prep with provided dispensing tip. Use provided brush to apply Vertise Flow to the occlusal surfaces with moderate pressure for 15-20 seconds to obtain a thin layer (<0.5mm). Remove excess material using the brush.
4. Light cure all coated surfaces for 20 seconds†. For A3.5 and Universal Opaque shades, light cure each increment for 40 seconds.
5. If necessary, add additional Vertise™ Flow and light cure for 20 seconds†. For A3.5 and Universal Opaque shades, light cure each increment for 40 seconds.
6. Check and adjust occlusion where necessary.

For Porcelain Repair

Sandblast the prep OR roughen the prep with a diamond bur. Rinse with water for 10 seconds, air-dry for 5 seconds.

1. Dispense Vertise™ Flow onto prep with provided dispensing tip. Use provided brush to apply Vertise™ Flow to the bonding surface(s) with moderate pressure for 15-20 seconds to obtain a thin layer (<0.5mm). Remove excess material around margins with the brush if necessary.
2. Light cure for 20 seconds.
3. If needed, an opaque can be applied and light-cured to hide the metallic color if metal is involved.

4. Build the restoration with more Vertise™ Flow or your composite of choice in increments of 2mm or less.
5. Light cure each increment of Vertise™ Flow for 20 seconds.† For A3.5 and Universal Opaque shades, light cure each increment for 40 seconds. For composites other than Vertise Flow, follow manufacturing instructions.
6. Finish and polish using a polishing system such as KerrHawe OptiDisc, Occlubrush, or HiLuster.

Note: Vertise™ Flow bendable brush is intended for single patient use only to prevent cross-contamination between patients.

† Use a curing light with output equal to or greater than 800 MW/cm² such as Kerr's Demi, L.E.Demetron II, or Optilux 501.

CAUTION: Uncured methacrylate resin may cause contact dermatitis and damage the pulp. Avoid contact with skin, eyes, and soft tissue. Wash thoroughly with water after contact.

Limited Warranty

Limitation of Kerr's Liability

Kerr's technical advice, whether verbal or in writing, is designed to assist dentists in using Kerr's product. The dentist assumes all risk and liability for damages arising out of the improper use of Kerr's product. In the event of a defect in material or workmanship, Kerr's liability is limited, at Kerr's option, to replacement of the defective product or part thereof, or reimbursement of the actual cost of the defective product. In order to take advantage of this limited warranty, the defective product must be returned to Kerr. In no event shall Kerr be liable for any indirect, incidental, or consequential damages. EXCEPT AS EXPRESSLY PROVIDED ABOVE, THERE ARE NO WARRANTIES, BY KERR, EXPRESS OR IMPLIED, INCLUDING WARRANTIES WITH RESPECT TO DESCRIPTION, QUALITY, OR FITNESS FOR A PARTICULAR PURPOSE.

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Limited Warranty

Limitation of Kerr's Liability

Kerr's technical advice, whether verbal or in writing, is designed to assist dentists in using Kerr's product. The dentist assumes all risk and liability for damages arising out of the improper use of Kerr's product. In the event of a defect in material or workmanship, Kerr's liability is limited, at Kerr's option, to replacement of the defective product or part thereof, or reimbursement of the actual cost of the defective product. In order to take advantage of this limited warranty, the defective product must be returned to Kerr. In no event shall Kerr be liable for any indirect, incidental, or consequential damages.

EXCEPT AS EXPRESSLY PROVIDED ABOVE, THERE ARE NO WARRANTIES, BY KERR, EXPRESS OR IMPLIED, INCLUDING WARRANTIES WITH RESPECT TO DESCRIPTION, QUALITY, OR FITNESS FOR A PARTICULAR PURPOSE.

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Kerr Corporation
Kerr Australia Pty. Ltd.
Sybron Dental Specialties Japan, Inc.